

Chemical Agent	CAS Number	Breakthrough Time	Protection Index	EN Std.	Materiaal	Naam	Brand
Acetic Acid, Glacial	64-19-7	2	0	374-3:1994	Latex	Conform 69-140	Ansell
Carbon disulfide	75-15-0	0	0	374-3:1994	Latex	Conform 69-140	Ansell
Formaldehyde, 35 %	50-00-0	325	5	374-3:1994	Latex	Conform 69-140	Ansell
Hydrofluoric Acid, 48 %	7664-39-3	33	2	374-3:1994	Latex	Conform 69-140	Ansell
Hydrogen Peroxide, 30%	7722-84-1	> 480	6	374-3:1994	Latex	Conform 69-140	Ansell
Nitric Acid, 50 %	7697-37-2	109	3	374-3:1994	Latex	Conform 69-140	Ansell
Sodium Hydroxide, 50 %	1310-73-2	> 480	6	374-3:1994	Latex	Conform 69-140	Ansell
Acetic Acid, Glacial	64-19-7	2	0	374-3:1994	Latex	Conform 69-150	Ansell
Acetone	67-64-1	< 5	0	374-3:1994	Latex	Conform 69-150	Ansell
Carbon disulfide	75-15-0	0	0	374-3:1994	Latex	Conform 69-150	Ansell
Dimethyl Sulfoxide	67-68-5	20	1	374-3:1994	Latex	Conform 69-150	Ansell
Formaldehyde, 35 %	50-00-0	396	5	374-3:1994	Latex	Conform 69-150	Ansell
Hydrofluoric Acid, 48 %	7664-39-3	24	2	374-3:1994	Latex	Conform 69-150	Ansell
Hydrogen Peroxide, 30%	7722-84-1	> 480	6	374-3:1994	Latex	Conform 69-150	Ansell
Nitric Acid, 50 %	7697-37-2	186	4	374-3:1994	Latex	Conform 69-150	Ansell
Sodium Hydroxide, 50 %	1310-73-2	> 480	6	374-3:1994	Latex	Conform 69-150	Ansell
Sulphuric acid, 50 %	7664-93-9	> 480	6	374-3:1994	Latex	Conform 69-150	Ansell
Acetic Acid, Glacial	64-19-7	5	0	374-3:1994	Latex	Gammex	Ansell
Dimethyl Sulfoxide	67-68-5	28	1	374-3:1994	Latex	Gammex	Ansell
Gamma Butyrolactone	96-48-0	10	1	374-3:1994	Latex	Gammex	Ansell
Methanol	67-56-1	3	0	374-3:1994	Latex	Gammex	Ansell
Propanol	71-23-8	2	0	374-3:1994	Latex	Gammex	Ansell
1-Methoxy-2-Propanol	107-98-2	7	0	374-3:1994	Nitril	TNT 92-500	Ansell
Acetic Acid, Glacial	64-19-7	1	0	374-3:1994	Nitril	TNT 92-500	Ansell
Ammonia, 25 %	7664-41-7	21	1	374-3:1994	Nitril	TNT 92-500	Ansell
Ammonium Hydroxide, 25 %	1336-21-6	67	3	374-3:1994	Nitril	TNT 92-500	Ansell
Ethanol	64-17-5	6	0	374-3:1994	Nitril	TNT 92-500	Ansell
Formaldehyde, 35 %	50-00-0	> 480	6	374-3:2003	Nitril	TNT 92-500	Ansell
Hydrogen Peroxide, 30%	7722-84-1	38	2	374-3:2003	Nitril	TNT 92-500	Ansell
Isopropanol	67-63-0	24	1	374-3:1994	Nitril	TNT 92-500	Ansell
Nitric Acid, 70 %	7697-37-2	< 5	0	374-3:1994	Nitril	TNT 92-500	Ansell
Perchloroethylene	127-18-4	< 5	0	374-3:1994	Nitril	TNT 92-500	Ansell
Sulphuric acid, 50 %	7664-93-9	> 480	6	374-3:1994	Nitril	TNT 92-500	Ansell
Triethylamine	121-44-8	9	0	374-3:1994	Nitril	TNT 92-500	Ansell
Xylene	1330-20-7	< 5	0	374-3:1994	Nitril	TNT 92-500	Ansell
1,1,1-trichloro-2-methyl-2-propyl alcohol	-	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
1,2-dichloroethane	107-06-2	20sec.	0	374-3:1994	Nitril	TNT 92-600	Ansell
1-Methoxy-2-Propanol	107-98-2	14.8	1	374-3:1994	Nitril	TNT 92-600	Ansell
Acetic Acid, Glacial	64-19-7	7	0	374-3:1994	Nitril	TNT 92-600	Ansell
Acetonitrile 73 % + Methyl Alcohol 25 %	-	1	0	374-3:1994	Nitril	TNT 92-600	Ansell
Acrylic Acid	79-10-7	71	3	374-3:1994	Nitril	TNT 92-600	Ansell
Allylchloride	107-05-1	75	3	374-3:1994	Nitril	TNT 92-600	Ansell
Ammonia, 25 %	7664-41-7	53	2	374-3:1994	Nitril	TNT 92-600	Ansell
Ammonium Hydroxide, 25 %	1336-21-6	101	3	374-3:1994	Nitril	TNT 92-600	Ansell
Benzyl Alcohol	100-51-6	11	1	374-3:1994	Nitril	TNT 92-600	Ansell
Bromochloromethane	74-97-5	104	3	374-3:1994	Nitril	TNT 92-600	Ansell
Butyl Alcohol	71-36-3	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Caffeine 1.6 %	58-08-2	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Carbon disulfide	75-15-0	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Chlorobutane	25154-42-1	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Chloroform	67-66-3	1	0	374-3:1994	Nitril	TNT 92-600	Ansell
Cyclohexane	110-82-7	> 480	6	374-3:2003	Nitril	TNT 92-600	Ansell
Cyclohexanone	108-94-1	8.2	0	374-3:1994	Nitril	TNT 92-600	Ansell
Dibromoethane	106-93-4	20sec.	0	374-3:1994	Nitril	TNT 92-600	Ansell
Dibromomethane	74-95-3	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Diesel fuel	68334-30-5	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Diethylamine	109-89-7	1	0	374-3:1994	Nitril	TNT 92-600	Ansell
Diethyl ether	60-29-7	< 40sec.	0	374-3:1994	Nitril	TNT 92-600	Ansell
Dimethylformamide	68-12-2	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Dimethyl Sulfoxide	67-68-5	5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Ethanol	64-17-5	8	0	374-3:1994	Nitril	TNT 92-600	Ansell
Ethanol, 70%	64-17-5	27	1	374-3:2003	Nitril	TNT 92-600	Ansell
Ethidium bromide in water (saturated)	1239-45-8	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Ethyl Acetate	141-78-6	< 40sec.	0	374-3:1994	Nitril	TNT 92-600	Ansell
Ethyl acetate 86 % + Methyl Alcohol 14 %	-	1	0	374-3:1994	Nitril	TNT 92-600	Ansell
Formaldehyde, 24.5 %	50-00-0	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Formaldehyde, 35 %	50-00-0	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Formaldehyde 4 % in Phosphatebuffer	50-00-0	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Gasoline	8006-61-9	98	3	374-3:1994	Nitril	TNT 92-600	Ansell
Glutaraldehyde, 50 %	111-30-8	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Heptane 98 % + 1-butyl alcohol 2 %	142-82-5	10	1	374-3:1994	Nitril	TNT 92-600	Ansell

Chemical Agent	CAS Number	Breakthrough Time	Protection Index	EN Std.	Materiaal	Naam	Brand
Heptane 98 % + 3-methyl-1-butyl alcohol	142-82-5	19	1	374-3:1994	Nitril	TNT 92-600	Ansell
Hexane	110-54-3	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Hydrochloric Acid, 37 %	7647-01-0	65	3	374-3:1994	Nitril	TNT 92-600	Ansell
Hydrofluoric Acid, 48 %	7664-39-3	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Hydrogen Bromide, 49 %	10035-10-6	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Hydrogen Peroxide, 30%	7722-84-1	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Iso-Octane	540-84-1	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Isopropanol	67-63-0	161	4	374-3:1994	Nitril	TNT 92-600	Ansell
Kerosene	64742-81-0	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Methanol	67-56-1	1	0	374-3:1994	Nitril	TNT 92-600	Ansell
Methyl ethyl ketone	78-93-3	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Methyl Isobutyl Ketone	108-10-1	1	0	374-3:1994	Nitril	TNT 92-600	Ansell
Methyl Sulfoxide 5 % in Citratebuffer	-	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Methylviolet 1 %	8004-87-3	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Nicotine	54-11-5	58	2	374-3:1994	Nitril	TNT 92-600	Ansell
Nitric Acid, 70 %	7697-37-2	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
n-Undecane	1120-21-4	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Peracetic acid, 39 %	79-21-0	9.8	0	374-3:1994	Nitril	TNT 92-600	Ansell
Perchloroethylene	127-18-4	14.9	1	374-3:1994	Nitril	TNT 92-600	Ansell
Potassium permanganate 5 %	7722-64-7	120	4	374-3:1994	Nitril	TNT 92-600	Ansell
Propanol	71-23-8	21	1	374-3:1994	Nitril	TNT 92-600	Ansell
Salicylic Acid 2 % in Peanut oil	-	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Sodium Hydroxide, 50 %	1310-73-2	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Sulphuric acid, 50 %	7664-93-9	> 480	6	374-3:1994	Nitril	TNT 92-600	Ansell
Sulphuric acid, 98 %	7664-93-9	2	0	374-3:1994	Nitril	TNT 92-600	Ansell
Tetrahydrofuran	109-99-9	< 5	0	374-3:1994	Nitril	TNT 92-600	Ansell
Toluene	108-88-3	1	0	374-3:1994	Nitril	TNT 92-600	Ansell
Triethylamine	121-44-8	41	2	374-3:1994	Nitril	TNT 92-600	Ansell
Xylene	1330-20-7	7.3	0	374-3:1994	Nitril	TNT 92-600	Ansell
1-Methoxy-2-Propanol	107-98-2	7	0	374-3:1994	Nitril	TNT Blue	Ansell
Acetic Acid, Glacial	64-19-7	7	0	374-3:1994	Nitril	TNT Blue	Ansell
Ammonia, 25 %	7664-41-7	17	1	374-3:1994	Nitril	TNT Blue	Ansell
Ammonium Hydroxide, 25 %	1336-21-6	52	2	374-3:1994	Nitril	TNT Blue	Ansell
Cyclohexane	110-82-7	18	1	374-3:2003	Nitril	TNT Blue	Ansell
Ethanol	64-17-5	6	0	374-3:1994	Nitril	TNT Blue	Ansell
Ethanol, 70%	64-17-5	12	1	374-3:2003	Nitril	TNT Blue	Ansell
Formaldehyde, 35 %	50-00-0	136	4	374-3:1994	Nitril	TNT Blue	Ansell
Hydrogen Peroxide, 30%	7722-84-1	> 480	6	374-3:1994	Nitril	TNT Blue	Ansell
Isopropanol	67-63-0	24	1	374-3:1994	Nitril	TNT Blue	Ansell
Nitric Acid, 70 %	7697-37-2	< 5	0	374-3:1994	Nitril	TNT Blue	Ansell
Perchloroethylene	127-18-4	< 5	0	374-3:1994	Nitril	TNT Blue	Ansell
Sulphuric acid, 50 %	7664-93-9	> 480	6	374-3:1994	Nitril	TNT Blue	Ansell
Triethylamine	121-44-8	15.9	1	374-3:1994	Nitril	TNT Blue	Ansell
Xylene	1330-20-7	< 5	0	374-3:1994	Nitril	TNT Blue	Ansell
2-Hydroxyethyl methacrylate	868-77-9	37	2	374-3:1994	Neopreen	Neotouch	Ansell
Acetone	67-64-1	0.3	0	374-3:1994	Neopreen	Neotouch	Ansell
Ammonia, 25 %	7664-41-7	9	0	374-3:1994	Neopreen	Neotouch	Ansell
Cefuroxim Sodium salt 15 g/l	-	> 480	6	374-3:1994	Neopreen	Neotouch	Ansell
Cyclohexane	110-82-7	< 5	0	374-3:2003	Neopreen	Neotouch	Ansell
Diethyl ether	60-29-7	0.3	0	374-3:1994	Neopreen	Neotouch	Ansell
Dimethylacetamide	127-19-5	4.9	0	374-3:1994	Neopreen	Neotouch	Ansell
Dimethylformamide	68-12-2	5	0	374-3:1994	Neopreen	Neotouch	Ansell
Dimethyl Sulfoxide	67-68-5	13	1	374-3:1994	Neopreen	Neotouch	Ansell
Ethanol, 70%	64-17-5	14	1	374-3:2003	Neopreen	Neotouch	Ansell
Ethanol, 96%	64-17-5	6.6	0	374-3:1994	Neopreen	Neotouch	Ansell
Formaldehyde, 35 %	50-00-0	> 480	6	374-3:2003	Neopreen	Neotouch	Ansell
Formaldehyde, 4%	50-00-0	> 480	6	374-3:2003	Neopreen	Neotouch	Ansell
Hydrochloric Acid, 37 %	7647-01-0	107	3	374-3:1994	Neopreen	Neotouch	Ansell
Hydrogen Peroxide, 30%	7722-84-1	> 480	6	374-3:1994	Neopreen	Neotouch	Ansell
Isopropanol	67-63-0	36	2	374-3:1994	Neopreen	Neotouch	Ansell
Methanol	67-56-1	9	0	374-3:1994	Neopreen	Neotouch	Ansell
Methyl ethyl ketone	78-93-3	0.3	0	374-3:1994	Neopreen	Neotouch	Ansell
Metronidazol solution 5 g/l	-	> 480	6	374-3:1994	Neopreen	Neotouch	Ansell
Nitric Acid, 70 %	7697-37-2	32	2	374-3:1994	Neopreen	Neotouch	Ansell
Perchloroethylene	127-18-4	< 5	0	374-3:1994	Neopreen	Neotouch	Ansell
Pyridine	110-86-1	0.7	0	374-3:1994	Neopreen	Neotouch	Ansell
Sodium Hydroxide, 50 %	1310-73-2	> 480	6	374-3:1994	Neopreen	Neotouch	Ansell
Sulphuric acid, 98 %	7664-93-9	7	0	374-3:1994	Neopreen	Neotouch	Ansell
Tetrahydrofuran	109-99-9	0.3	0	374-3:1994	Neopreen	Neotouch	Ansell
Toluene	108-88-3	0.7	0	374-3:1994	Neopreen	Neotouch	Ansell
Triethylamine	121-44-8	< 5	0	374-3:1994	Neopreen	Neotouch	Ansell

Chemical Agent	CAS Number	Breakthrough Time	Protection Index	EN Std.	Materiaal	Naam	Brand
Xylene	1330-20-7	< 5	0	374-3:1994	Neopreen	Neotouch	Ansell
Acetone	67-64-1	2.2	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Acetonitrile	75-05-8	8.7	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Butyl Alcohol	71-36-3	> 480	6	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Cyclohexane	110-82-7	7	0	374-3:2003	Neopreen	Derma Shield 73-701	Ansell
Diethyl ether	60-29-7	0.9	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Dimethylacetamide	127-19-5	7.3	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Ethanol	64-17-5	78	3	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Ethanol, 70%	64-17-5	78	3	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Ethanol 95% + Dettol 5%	-	26	1	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Formaldehyde, 35 %	50-00-0	> 480	6	374-3:2003	Neopreen	Derma Shield 73-701	Ansell
Formaldehyde, 4%	50-00-0	> 480	6	374-3:2003	Neopreen	Derma Shield 73-701	Ansell
Hydrochloric Acid, 37 %	7647-01-0	> 480	6	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Isopropanol	67-63-0	> 480	6	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Methanol	67-56-1	38	2	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Methyl ethyl ketone	78-93-3	1.25	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Nitric Acid, 70 %	7697-37-2	> 480	6	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Pyridine	110-86-1	1.5	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Sodium Hydroxide, 50 %	1310-73-2	> 480	6	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Sulphuric acid, 95 %	7664-93-9	229	4	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Tetrahydrofuran	109-99-9	0.3	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Toluene	108-88-3	1	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Triethylamine	121-44-8	< 5	0	374-3:1994	Neopreen	Derma Shield 73-701	Ansell
Hydrogen Peroxide, 30%	7722-84-1	> 480	6	374-3:1994	Vinyl	DuraTouch	Ansell
Isopropanol	67-63-0	< 5	0	374-3:1994	Vinyl	DuraTouch	Ansell
Sodium Hydroxide, 50 %	1310-73-2	> 480	6	374-3:1994	Vinyl	DuraTouch	Ansell
Minimaal risico	-	nvt	-	Cat-1	Vinyl	Soft Star Vinyl	Samurai
Minimaal risico	-	nvt	-	Cat-1	Vinyl	Pro-Dispro Vinyl P	Artelli
Minimaal risico	-	nvt	-	Cat-1	Vinyl	Pro-Dispro Vinyl P BL	Artelli
Minimaal risico	-	nvt	-	Cat-1	Vinyl	Disp. Vinyl P blauw	GSA
Minimaal risico	-	nvt	-	Cat-1	Vinyl	Disp. Vinyl P wit	GSA
Minimaal risico	-	nvt	-	Cat-1	Vinyl	Intacto Vinyl	Intacto
Minimaal risico	-	nvt	-	Cat-1	Vinyl	Intacto Vinyl Blauw	Intacto
Dimethylformamide	68-12-2	> 480	6	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Ethyl acetate	141-78-6	> 480	6	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Hexane n-	110-54-3	> 60	3	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Methanol	67-56-1	> 480	6	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Sodium hydroxide 40%	1310-73-2	> 480	6	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Nitrobenzene	98-95-3	< 10	0	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Tetrahydrofurane	109-99-9	> 10	1	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Sulphuric acid 90%	7664-93-9	> 10	1	374-3:1994	Nitrile	Disp Nitrile BL	GSA
Dimethylformamide	68-12-2	> 480	6	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Ethyl acetate	141-78-6	> 480	6	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Hexane n-	110-54-3	> 60	3	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Methanol	67-56-1	> 480	6	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Sodium hydroxide 40%	1310-73-2	> 480	6	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Nitrobenzene	98-95-3	< 10	0	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Tetrahydrofurane	109-99-9	> 10	1	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Sulphuric acid 90%	7664-93-9	> 10	1	374-3:1994	Nitrile	Disp Nitrile P BL	GSA
Dimethylformamide	68-12-2	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Ethyl acetate	141-78-6	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Hexane n-	110-54-3	> 60	3	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Methanol	67-56-1	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Sodium hydroxide 40%	1310-73-2	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Nitrobenzene	98-95-3	< 10	0	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Tetrahydrofurane	109-99-9	> 10	1	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Sulphuric acid 90%	7664-93-9	> 10	1	374-3:1994	Nitrile	Pro-Dispro Nitril P	Artelli
Dimethylformamide	68-12-2	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Ethyl acetate	141-78-6	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Hexane n-	110-54-3	> 60	3	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Methanol	67-56-1	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Sodium hydroxide 40%	1310-73-2	> 480	6	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Nitrobenzene	98-95-3	< 10	0	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Tetrahydrofurane	109-99-9	> 10	1	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Sulphuric acid 90%	7664-93-9	> 10	1	374-3:1994	Nitrile	Pro-Dispro Nitril PF	Artelli
Dimethylformamide	68-12-2	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile	Samurai
Ethyl acetate	141-78-6	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile	Samurai
Hexane n-	110-54-3	> 60	3	374-3:1994	Nitrile	Soft Star Nitrile	Samurai
Methanol	67-56-1	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile	Samurai
Sodium hydroxide 40%	1310-73-2	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile	Samurai
Nitrobenzene	98-95-3	< 10	0	374-3:1994	Nitrile	Soft Star Nitrile	Samurai

Chemical Agent	CAS Number	Breakthrough Time	Protection Index	EN Std.	Materiaal	Naam	Brand
Tetrahydrofurane	109-99-9	> 10	1	374-3:1994	Nitrile	Soft Star Nitrile	Samurai
Sulphuric acid 90%	7664-93-9	> 10	1	374-3:1994	Nitrile	Soft Star Nitrile	Samurai
Dimethylformamide	68-12-2	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai
Ethyl acetate	141-78-6	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai
Hexane n-	110-54-3	> 60	3	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai
Methanol	67-56-1	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai
Sodium hydroxide 40%	1310-73-2	> 480	6	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai
Nitrobenzene	98-95-3	< 10	0	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai
Tetrahydrofurane	109-99-9	> 10	1	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai
Sulphuric acid 90%	7664-93-9	> 10	1	374-3:1994	Nitrile	Soft Star Nitrile P	Samurai

Ansell Permeatie index	
BTT minuten	
Rood	< 10
Oranje	10 - 60
Geel	60 - 240
Groen	> 240

Permeatie index EN 374-3 BTT	
Minuten	Index
< 10	0
> 10	1
> 30	2
> 60	3
> 120	4
> 240	5
> 480	6

Geen geg.	wit
-----------	-----

	Chemical Agent	CAS Number	Breakthrough Time	Protection Index	EN Std.	Materiaal	Naam	Brand
--	----------------	------------	-------------------	------------------	---------	-----------	------	-------

